

## PRODUCT BROCHURE

---

# KERN 347-057

Date: 10.03.2026

ENGLISH



Individual weight, OIML Class M1, 20 g, Slot weight (Type B), stainless steel fine turned

---

### CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)



Scan here for more  
information

## DESCRIPTION

- Individual weight, OIML Class M1, 20 g, Slot weight (Type B), stainless steelfine turned (OIML)

## SPECIFICATIONS

### Article codes, product groups & Status information

|                       |                                |
|-----------------------|--------------------------------|
| Model article-code    | 347-057                        |
| Model-Serie           | 347-057                        |
| Assortment            | KERN                           |
| Product Type          | Test weights                   |
| Product group         | OIML M1 Test weights           |
| Product usage type    | Main item;Accessories/Services |
| EAN code (Model)      | 4045761366657                  |
| Customs tariff number | 84239090                       |
| Country of Origin     | DE                             |

### Technical data - Primary

|                  |        |
|------------------|--------|
| Nominal value    | 20 g   |
| Tolerance (OIML) | 2,5 mg |

### Technical data - Construction

|          |                             |
|----------|-----------------------------|
| Material | stainless steel fine turned |
|----------|-----------------------------|

### Technical data - Measuring system

|               |      |
|---------------|------|
| Nominal value | 20 g |
| OIML class    | M1   |

|                              |        |
|------------------------------|--------|
| OIML conform                 | yes    |
| Tolerance (OIML)             | 2,5 mg |
| Density (kg/m <sup>3</sup> ) | 8000   |
| Volume (cm <sup>3</sup> )    | 2,381  |

### Technical data - Packing & Shipping

|                                            |                |
|--------------------------------------------|----------------|
| Dimensions packaging (W×D×H)               | 38×38×4 mm     |
| Gross weight                               | 0.05 kg        |
| Net weight                                 | 0.05 kg        |
| Shipping method                            | Parcel service |
| Packaging component - by weight - plastics | 0 g            |
| Delivery time                              | 1 d            |

### Optional Services

|                   |         |
|-------------------|---------|
| DAkkS Calibration | 962-635 |
|-------------------|---------|

## FUNCTIONS

|          |
|----------|
| Standard |
|----------|



|        |
|--------|
| Option |
|--------|



## ACCESSORIES

| Model                              | Description                     |
|------------------------------------|---------------------------------|
| <a href="#"><u>347-137-100</u></a> | Supporting Bar KERN 347-137-100 |
| <a href="#"><u>347-147-100</u></a> | Supporting Bar KERN 347-147-100 |
| <a href="#"><u>347-117-100</u></a> | Supporting Bar KERN 347-117-100 |
| <a href="#"><u>315-242</u></a>     | Tweezers KERN 315-242           |

## SERVICES

| Model                           | Description                       |
|---------------------------------|-----------------------------------|
| <a href="#"><u>962-635</u></a>  | DAkkS Calibration KERN 962-635    |
| <a href="#"><u>962-635R</u></a> | DAkkS Recalibration KERN 962-635R |