

PRODUCT BROCHURE

KERN 317-02

Date: 12.03.2026

ENGLISH



Individual weight, OIML Class E2, 2 g, Knob, stainless steel polished

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

- Individual weight, OIML Class E2, 2 g, Knob, stainless steel polished (OIML)
- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section
- For the verification you need a box. Please add a box for single weights.

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	317-02
Model-Serie	317
Assortment	KERN
Product Type	Test weights
Product group	OIML E2 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	404576 1312500
Predecessor Model	317-02-2020a; 317-02-2022i
Customs tariff number	84239090
Country of Origin	CN

Technical data - Primary

Nominal value	2 g
Tolerance (OIML)	0,04 mg
Case aluminium	317-020-600
Case plastic	317-020-400
Case wood	317-020-100

Technical data - Construction

Material	stainless steel polished
----------	--------------------------

Technical data - Measuring system

Nominal value	2 g
OIML class	E2
OIML conform	yes
Tolerance (OIML)	0,04 mg
Density (kg/m ³)	7950
Volume (cm ³)	0,25

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	6×6×10 mm
Gross weight	0.05 kg
Net weight	0.05 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible yes

Optional Services

DAkkS Calibration 962-332

FUNCTIONS

Standard



Option



ACCESSORIES

Model	Description
<u>317-020-100</u>	Wooden Box for Set of Weights KERN 317-020-100
<u>317-020-400</u>	Plastic Weight Box KERN 317-020-400
<u>317-020-600</u>	Aluminium Protective Box KERN 317-020-600
<u>315-242</u>	Tweezers KERN 315-242

SERVICES

Model	Description
<u>962-332</u>	DAkkS Calibration KERN 962-332
<u>962-332R</u>	DAkkS Recalibration KERN 962-332R
<u>921-332</u>	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-332

RELATED ARTICLES DATA

Model name	Nominal value	Tolerance (OIML)	Case aluminium	Case plastic	Case wood
<u>317-09</u>	500 g	0,8 mg	317-090-600	317-090-400	317-090-100
<u>317-11</u>	1 kg	1,6 mg	317-110-600	317-110-400	317-110-100

Model name	Nominal value	Tolerance (OIML)	Case aluminium	Case plastic	Case wood
<u>317-12</u>	2 kg	3 mg	317-120-600	317-120-400	317-120-100
<u>317-13</u>	5 kg	8 mg	317-130-600	317-130-400	317-130-100
<u>317-14</u>	10 kg	16 mg	317-140-600	317-140-400	317-140-100
<u>317-15</u>	20 kg	30 mg	317-150-600	-	317-150-100
<u>317-16</u>	50 kg	80 mg	317-160-600	-	317-160-100
<u>317-01</u>	1 g	0,03 mg	317-010-600	317-020-400	317-010-100
<u>317-02</u>	2 g	0,04 mg	317-020-600	317-020-400	317-020-100
<u>317-03</u>	5 g	0,05 mg	317-030-600	317-030-400	317-030-100
<u>317-04</u>	10 g	0,06 mg	317-040-600	317-040-400	317-040-100
<u>317-05</u>	20 g	0,08 mg	317-050-600	317-050-400	317-050-100
<u>317-06</u>	50 g	0,1 mg	317-060-600	317-060-400	317-060-100
<u>317-07</u>	100 g	0,16 mg	317-070-600	317-070-400	317-070-100
<u>317-08</u>	200 g	0,3 mg	317-080-600	317-080-400	317-080-100