


PRODUCT BROCHURE

KERN 347-087-N00

Date: 23.04.2026

ENGLISH



Individual weight (Newton), 2 N, Slot weight (Type B), stainless steel fine turned

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	347-087-N00
Model-Serie	347-7-N00
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	404576 1366725
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Tolerance (OIML)	10 mg
------------------	-------

Technical data - Construction

Material	stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

OIML class	M1
Density (kg/m ³)	8000
Volume (cm ³)	23,8095

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	60×60×14.1 mm
------------------------------	---------------

Gross weight	0.20 kg
Net weight	0.20 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	6 d

Optional Services

DAkkS Calibration	962-638N
-------------------	----------

FUNCTIONS

Standard



Option



ACCESSORIES

Model	Description
347-147-100-N00	Supporting Bar KERN 347-147-100-N00
347-137-100-N00	Supporting Bar KERN 347-137-100-N00
347-117-100-N00	Supporting Bar KERN 347-117-100-N00
315-242	Tweezers KERN 315-242

SERVICES

Model	Description
962-638RN	DAkkS Recalibration KERN 962-638RN
962-638N	DAkkS Calibration KERN 962-638N

RELATED ARTICLES DATA

Model name	Tolerance (OIML)
347-037-N00	1,6 mg
347-047-N00	2 mg
347-057-N00	2,5 mg
347-067-N00	3 mg
347-147-N00	500 mg
347-147-100-N00	500 mg
347-137-N00	250 mg
347-137-100-N00	250 mg
347-127-N00	100 mg
347-117-N00	50 mg
347-117-100-N00	50 mg
347-097-N00	25 mg
347-087-N00	10 mg
347-077-N00	5 mg