


KERN 346-055-009

Weight container for rectangular weights, calibrated to OIML class M1



Article codes, product groups & Status information

Model article-code	346-055-009
Model-Serie	346
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main article
EAN code (Model)	4045761225091
Planned start date	2019-1HJ
Customs tariff number	84239090

Net weight	50 kg
Shipping method	Freight
Packaging component - by weight - plastics	0 g
Delivery time	2 d

Optional Services

DAkkS Calibration	962-646
-------------------	---------

Technical data - Primary

Nominal value	50 kg
Tolerance (OIML)	2,5 g

Technical data - Construction

Dimensions (W×D×H)	1000×500×500 mm
Material	Stainless steel

Technical data - Measuring system

Nominal value	50 kg
OIML class	M1
OIML conform	Yes
Tolerance (OIML)	2,5 mg

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	1040×635×370 mm
Gross weight	50 kg

Functions

KERN 346-055-009

Weight container for rectangular weights, calibrated to OIML class M1



Standard

