

PRODUCT BROCHURE

KERN 314-086

Date: 10.03.2026

ENGLISH



Weight set, E2, 1 g - 5 kg, Knob, stainless steel polished, in aluminium case

CONTACT


KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	314-086
Model-Serie	314
Assortment	KERN
Product Type	Test weights
Product group	OIML E2 Test weights
Product usage type	Main article
EAN code (Model)	4045761223165
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 g - 5 kg
Case aluminium	313-080-600
Weight set case	in aluminium case

Technical data - Construction

Material	stainless steel polished
Is weight set	yes
Type of case	in aluminium case

Technical data - Measuring system

Nominal value	1 g - 5 kg
OIML class	E2
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	460×260×200 mm
Gross weight	18 kg
Net weight	13 kg
Shipping method	Parcel service

Packaging component - 5.4 kg
by weight - plastics

Delivery time 1 d

Technical data - Verification

Verification possible yes

Optional Services

DAkkS Calibration 962-321

FUNCTIONS

Standard



Option



SERVICES

Model	Description
<u>962-321</u>	DAkkS Calibration KERN 962-321
<u>962-321R</u>	DAkkS Recalibration KERN 962-321R
<u>921-321</u>	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-321
<u>920-716</u>	Marking (Last digits of serial number) KERN 920-716

SPARE PARTS

Model	Description
<u>313-080-600</u>	Aluminium Protective Box KERN 313-080-600
<u>315-242</u>	Tweezers KERN 315-242