

KERN 357-14

Individual weight, OIML Class M2, 10 kg, Knob, stainless steel fine turned



Article codes, product groups & Status information

Model article-code	357-14
Model-Serie	357-14
Assortment	KERN
Product Type	Test weights
Product group	OIML M2 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761289062
Planned start date	2022-1HJ
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	10 kg
Tolerance (OIML)	1,6 g
Case aluminium	317-140-600
Case plastic	347-140-400
Case wood	337-140-200

Technical data - Construction

Material	Stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	10 kg
OIML class	M2
OIML conform	Yes
Tolerance (OIML)	1,6 g

Functions

Standard

Density (kg/m ³)	7.950
Volume (cm ³)	1.190,4762

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	115×115×198 mm
Gross weight	10 kg
Net weight	10 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-644
-------------------	---------

KERN 357-14

Individual weight, OIML Class M2, 10 kg, Knob, stainless steel fine turned



Option

