

## PRODUCT BROCHURE

---

# KERN 366-77

Date: 29.04.2025

ENGLISH



Individual weight, OIML Class M3, 10 kg, Block Eco, cast iron lacquered

---

### CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)



Scan here for more  
information

## SPECIFICATIONS

### Article codes, product groups & Status information

|                       |                                |
|-----------------------|--------------------------------|
| Model article-code    | 366-77                         |
| Model-Serie           | 366                            |
| Assortment            | KERN                           |
| Product Type          | Test weights                   |
| Product group         | OIML M3 Test weights           |
| Product usage type    | Main item;Accessories/Services |
| EAN code (Model)      | 404576 124 1831                |
| Customs tariff number | 84239090                       |

### Technical data - Primary

|                  |             |
|------------------|-------------|
| Nominal value    | 10 kg       |
| Tolerance (OIML) | 5 g         |
| Case aluminium   | 346-070-600 |

### Technical data - Construction

|                    |                     |
|--------------------|---------------------|
| Dimensions (W×D×H) | 190×95×109 mm       |
| Material           | cast iron lacquered |

### Technical data - Measuring system

|                              |           |
|------------------------------|-----------|
| Nominal value                | 10 kg     |
| OIML class                   | M3        |
| OIML conform                 | yes       |
| Tolerance (OIML)             | 5 g       |
| Density (kg/m <sup>3</sup> ) | 7100      |
| Volume (cm <sup>3</sup> )    | 1408,4507 |

### Technical data - Packing & Shipping

|                                            |                |
|--------------------------------------------|----------------|
| Dimensions packaging (W×D×H)               | 190×95×109 mm  |
| Gross weight                               | 10 kg          |
| Net weight                                 | 10 kg          |
| Shipping method                            | Parcel service |
| Packaging component - by weight - plastics | 0 g            |
| Delivery time                              | 1 d            |

### Technical data - Verification

|                       |     |
|-----------------------|-----|
| Verification possible | yes |
|-----------------------|-----|

### Optional Services

|                   |         |
|-------------------|---------|
| DAkkS Calibration | 962-644 |
|-------------------|---------|

## FUNCTIONS

|          |
|----------|
| Standard |
|----------|



## ACCESSORIES

| Model       | Description                               |
|-------------|-------------------------------------------|
| 317-280     | Gloves KERN 317-280                       |
| 346-070-600 | Aluminium Protective Box KERN 346-070-600 |

## SERVICES

| Model    | Description                                                            |
|----------|------------------------------------------------------------------------|
| 962-644  | DAkkS Calibration KERN 962-644                                         |
| 962-644R | DAkkS Recalibration KERN 962-644R                                      |
| 952-644  | Verification KERN 952-644                                              |
| 920-30   | Marking with stroke number KERN 920-30                                 |
| 920-301  | Marking with stroke number (Last digits of serial number) KERN 920-301 |

## RELATED ARTICLES DATA

| Model name   | Nominal value | Tolerance (OIML) | Case aluminium |
|--------------|---------------|------------------|----------------|
| 366-94       | 1 kg          | 500 mg           | -              |
| 366-95       | 2 kg          | 1 g              | -              |
| 366-96       | 5 kg          | 2,5 g            | -              |
| 366-97       | 10 kg         | 5 g              | -              |
| 366-86       | 5 kg          | 2,5 g            | 346-060-600    |
| 366-87       | 10 kg         | 5 g              | 346-070-600    |
| 366-88       | 20 kg         | 10 g             | 346-080-600    |
| 366-89       | 50 kg         | 25 g             | -              |
| 366-91       | 100 g         | 50 mg            | -              |
| 366-92       | 200 g         | 100 mg           | -              |
| 366-93       | 500 g         | 250 mg           | -              |
| 366-76       | 5 kg          | 2,5 g            | 346-060-600    |
| 366-77       | 10 kg         | 5 g              | 346-070-600    |
| 366-78       | 20 kg         | 10 g             | 346-080-600    |
| 366-79       | 50 kg         | 25 g             | -              |
| 366-76-2025e | 5 kg          | 2,5 g            | -              |