

PRODUCT BROCHURE

KERN 334-07

Date: 09.03.2026

ENGLISH



Weight set, F2, 1 g - 2 kg, Knob, stainless steel fine turned, in wooden case

CONTACT


KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	334-07
Model-Serie	334
Assortment	KERN
Product Type	Test weights
Product group	OIML F2 Test weights
Product usage type	Main article
EAN code (Model)	4045761005976
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 g - 2 kg
Case wood	353-470-200
Weight set case	in wooden case

Technical data - Construction

Material	stainless steel fine turned
Is weight set	yes
Number of weights in set	15
Type of case	in wooden case

Technical data - Measuring system

Nominal value	1 g - 2 kg
OIML class	F2
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	305×180×140 mm
Gross weight	10 kg
Net weight	10 kg

Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-420
-------------------	---------

FUNCTIONS

Standard



Option



SERVICES

Model	Description
<u>962-420</u>	DAkkS Calibration KERN 962-420
<u>962-420R</u>	DAkkS Recalibration KERN 962-420R
<u>921-420</u>	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-420
<u>920-715</u>	Marking (Last digits of serial number) KERN 920-715

SPARE PARTS

Model	Description
<u>353-470-200</u>	Box of Weights KERN 353-470-200
<u>315-242</u>	Tweezers KERN 315-242