

KERN 328-39

Milligram weights, F1, 500 mg, Wire, stainless steel polished



Article codes, product groups & Status information

Model article-code	328-39
Model-Serie	328
Assortment	KERN
Product Type	Test weights
Product group	OIML F1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761478664
Predecessor Model	328-39-2.024a
Planned start date	2025-2HJ
Customs tariff number	84239090
Country of Origin	CN

Technical data - Primary

Nominal value	500 mg
Tolerance (OIML)	0,08 mg
Case aluminium	317-009-600
Case plastic	347-009-400
Case wood	338-090-200

Technical data - Construction

Dimensions (W×D×H)	15×15×4 mm
Material	Stainless steel polished

Technical data - Measuring system

Nominal value	500 mg
OIML class	F1

Functions

Standard

OIML conform	Yes
Tolerance (OIML)	0,08 mg
Density (kg/m ³)	8.000
Volume (cm ³)	0,0625

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	15×15×4 mm
Gross weight	0,05 kg
Net weight	0,05 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-459
-------------------	---------

KERN 328-39

Milligram weights, F1, 500 mg, Wire, stainless steel polished



Option

