

PRODUCT BROCHURE

OPTICS OZB-A4532

Date: 30.04.2025

ENGLISH



Double gooseneck 490 mm; in conjunction with OZB-A4512

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

External lighting unit OZB-A4532

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	OZB-A4532
Model-Serie	OZB-IF
Assortment	Optics
Product Type	Microscopes
Product group	External Illumination
Product usage type	Spare part;Accessories/Services
EAN code (Model)	4045761170834
Customs tariff number	90119000

Gross weight 1.00 kg

Net weight 0.70 kg

Shipping method Parcel service

Packaging component -
by weight - cardboard 150 g

Packaging component -
by weight - plastics 150 g

Delivery time 1 d

FUNCTIONS

Standard



ENGLISH

Technical data - Primary

Illumination type	Halogen
-------------------	---------

Technical data - Construction

Fibre length	490 mm
--------------	--------

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Stotage Temp Max	40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	280×185×50 mm
---------------------------------	---------------