

KERN 366-79

Individual weight, OIML Class M3, 50 kg, Block Eco, cast iron lacquered



Article codes, product groups & Status information

Model article-code	366-79
Model-Serie	366
Assortment	KERN
Product Type	Test weights
Product group	OIML M3 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761241855
Planned start date	2020-1HJ
Customs tariff number	84239090

Technical data - Primary

Nominal value	50 kg
Tolerance (OIML)	25 g

Technical data - Construction

Dimensions	310×155×192 mm
Material	Cast iron lacquered

Technical data - Measuring system

Nominal value	50 kg
OIML class	M3
OIML conform	Yes
Tolerance (OIML)	25 g
Density (kg/m ³)	7100
Volume (cm ³)	7042,2535

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	310×155×192 mm
Gross weight	50 kg
Net weight	50 kg
Shipping method	Freight
Packaging component - by weight - plastics	0 g
Delivery time	2 d

Technical data - Verification

Verification possible	Yes
-----------------------	-----

Optional Services

DAkKS Calibration	962-646
-------------------	---------

Functions

KERN 366-79

Individual weight, OIML Class M3, 50 kg, Block Eco, cast iron lacquered



Standard

