

## PRODUCT BROCHURE

---

# KERN 312-074

Date: 09.03.2026

ENGLISH



Weight set, E2, 1 g - 2 kg, Compact form, stainless steel polished, in plastic case

---

### CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)



Scan here for more  
information

## PRODUCT SHOWCASE



## SPECIFICATIONS

### Article codes, product groups & Status information

|                       |                      |
|-----------------------|----------------------|
| Model article-code    | 312-074              |
| Model-Serie           | 312                  |
| Assortment            | KERN                 |
| Product Type          | Test weights         |
| Product group         | OIML E2 Test weights |
| Product usage type    | Main article         |
| EAN code (Model)      | 4045761117037        |
| Customs tariff number | 84239090             |
| Country of Origin     | DE                   |

### Technical data - Primary

|                 |                 |
|-----------------|-----------------|
| Nominal value   | 1 g - 2 kg      |
| Case plastic    | 313-082-400     |
| Weight set case | in plastic case |

### Technical data - Construction

|                          |                          |
|--------------------------|--------------------------|
| Material                 | stainless steel polished |
| Is weight set            | yes                      |
| Number of weights in set | 15                       |
| Type of case             | in plastic case          |

### Technical data - Measuring system

|               |            |
|---------------|------------|
| Nominal value | 1 g - 2 kg |
| OIML class    | E2         |
| OIML conform  | yes        |

### Technical data - Packing & Shipping

|                              |                |
|------------------------------|----------------|
| Dimensions packaging (W×D×H) | 350×320×180 mm |
| Gross weight                 | 8 kg           |
| Net weight                   | 8 kg           |
| Shipping method              | Parcel service |

Packaging component - 0 g  
by weight - plastics

Delivery time 1 d

## Technical data - Verification

Verification possible no

## Optional Services

DAkkS Calibration 962-320

## FUNCTIONS

Standard



Option



## SERVICES

| Model                           | Description                                          |
|---------------------------------|------------------------------------------------------|
| <a href="#"><u>962-320</u></a>  | DAkkS Calibration KERN 962-320                       |
| <a href="#"><u>962-320R</u></a> | DAkkS Recalibration KERN 962-320R                    |
| <a href="#"><u>921-320</u></a>  | Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-320 |
| <a href="#"><u>920-715</u></a>  | Marking (Last digits of serial number) KERN 920-715  |

## SPARE PARTS

| Model                              | Description                         |
|------------------------------------|-------------------------------------|
| <a href="#"><u>313-082-400</u></a> | Plastic Weight Box KERN 313-082-400 |
| <a href="#"><u>315-242</u></a>     | Tweezers KERN 315-242               |