



The digital specialist for measurement, counting and archiving – free of charge with all KERN microscope cameras

Features

- As well as the streaming function for the object to be viewed, the software OXM 902 also offers you an image snapshot function, as well as a video function
- Various measuring functions such as, for example, functions for measuring distance, surfaces and angles and a manual counting function are available. In addition there are extensive image processing and documentation functions available, and of course an exporting function to Microsoft® Office applications Word and Excel
- With the display settings you can display different measurements, grid sizes, scales and rulers for the very best measuring procedure

- Automatic white balance and automatic contrast adjustment can be performed quickly and easily, which enables efficient working procedures
- The following highlights are included:
 - Image Stitching
 - Image Stacking
 - DShow and TWAIN support
 - Software development kit
- With this software it is possible to operate all available KERN microscope cameras

Technical data

- Can be used for Microsoft® Windows 7, 8, 10 and 11
- Depending on the language setting of your Windows operating system the software KERN microscope VIS Basic & Pro will be identified and installed in the current language. This can be changed manually at any time
- The software is available in the following languages: German, English, Spanish, Italian, French, Portuguese, Polish, Russian, Turkish, Chinese, Japanese, Korean
- Please refer to the documentation for the software in the download area of our website www.kern-sohn.com for further details