

PRODUCT BROCHURE

KERN 317-14

Date: 09.03.2026

ENGLISH



Individual weight, OIML Class E2, 10 kg, Knob, stainless steel polished

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

- Individual weight, OIML Class E2, 10 kg, Knob, stainless steel polished (OIML)
- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section
- For the verification you need a box. Please add a box for single weights.

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	317-14
Model-Serie	317
Assortment	KERN
Product Type	Test weights
Product group	OIML E2 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	404576 1003767
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	10 kg
Tolerance (OIML)	16 mg
Case aluminium	317-140-600
Case plastic	317-140-400
Case wood	317-140-100

Technical data - Construction

Material	stainless steel polished
----------	--------------------------

Technical data - Measuring system

Nominal value	10 kg
OIML class	E2
OIML conform	yes
Tolerance (OIML)	16 mg
Density (kg/m ³)	7950
Volume (cm ³)	1250

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	100×100×185 mm
Gross weight	10 kg
Net weight	10 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	yes
-----------------------	-----

Optional Services

DAkkS Calibration	962-344
-------------------	---------

FUNCTIONS

Standard



Option



ACCESSORIES

Model	Description
317-280	Gloves KERN 317-280
317-290	Gloves KERN 317-290
317-140-100	Wooden Box for Set of Weights KERN 317-140-100
317-140-400	Plastic Weight Box KERN 317-140-400
317-140-600	Aluminium Protective Box KERN 317-140-600

SERVICES

Model	Description
962-344	DAkkS Calibration KERN 962-344
962-344R	DAkkS Recalibration KERN 962-344R
921-344	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-344