

KERN 317-06

Individual weight, OIML Class E2, 50 g, Knob, stainless steel polished



Article codes, product groups & Status information

Model article-code	317-06
Model-Serie	317
Assortment	KERN
Product Type	Test weights
Product group	OIML E2 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761474949
Predecessor Model	317-6-2.020a; 317-6-2.022i
Planned start date	2022-2HJ
Customs tariff number	84239090
Country of Origin	CN

Technical data - Primary

Nominal value	50 g
Tolerance (OIML)	0,1 mg
Case aluminium	317-060-600
Case plastic	317-060-400
Case wood	317-060-100

Technical data - Construction

Material	Stainless steel polished
----------	--------------------------

Technical data - Measuring system

Nominal value	50 g
OIML class	E2
OIML conform	Yes

Functions

Standard

Tolerance (OIML)	0,1 mg
Density (kg/m ³)	7.950
Volume (cm ³)	6,25

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	18×18×28 mm
Gross weight	0,05 kg
Net weight	0,05 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkKS Calibration	962-336
-------------------	---------

KERN 317-06

Individual weight, OIML Class E2, 50 g, Knob, stainless steel polished



Option

