

PRODUCT BROCHURE

KERN 344-034

Date: 14.05.2026

ENGLISH



Weight set, M1, 1 g - 100 g, Knob, stainless steel fine turned, in plastic case

CONTACT


KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	344-034
Model-Serie	344
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main article
EAN code (Model)	404576 10079 18
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 g - 100 g
Case plastic	313-052-400
Weight set case	in plastic case

Technical data - Construction

Material	stainless steel fine turned
Is weight set	yes
Number of weights in set	9
Type of case	in plastic case

Technical data - Measuring system

Nominal value	1 g - 100 g
OIML class	M1
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	250×210×75 mm
Gross weight	0.60 kg
Net weight	0.60 kg

Shipping method Parcel service

Packaging component -
by weight - plastics 0 g

Delivery time 1 d

Optional Services

DAkkS Calibration 962-616

FUNCTIONS

Standard



Option

DAkkS

SERVICES

Model	Description
962-616	DAkkS Calibration KERN 962-616
962-616R	DAkkS Recalibration KERN 962-616R
920-711	Test weight marking KERN 920-711

SPARE PARTS

Model	Description
313-052-400	Plastic Weight Box KERN 313-052-400
315-242	Tweezers KERN 315-242

RELATED ARTICLES DATA

Model name	Nominal value	Case aluminium	Case plastic	Case wood	Weight set case
344-02	1 g - 50 g	-	-	353-420-200	in wooden case
344-024	1 g - 50 g	-	313-052-400	-	in plastic case

Model name	Nominal value	Case aluminium	Case plastic	Case wood	Weight set case
<u>344-03</u>	1 g - 100 g	-	-	353-430-200	in wooden case
<u>344-034</u>	1 g - 100 g	-	313-052-400	-	in plastic case
<u>344-04</u>	1 g - 200 g	-	-	353-440-200	in wooden case
<u>344-044</u>	1 g - 200 g	-	313-052-400	-	in plastic case
<u>344-05</u>	1 g - 500 g	-	-	353-450-200	in wooden case
<u>344-054</u>	1 g - 500 g	-	313-052-400	-	in plastic case
<u>344-06</u>	1 g - 1 kg	-	-	353-460-200	in wooden case
<u>344-064</u>	1 g - 1 kg	-	313-082-400	-	in plastic case
<u>344-07</u>	1 g - 2 kg	-	-	353-470-200	in wooden case
<u>344-074</u>	1 g - 2 kg	-	313-082-400	-	in plastic case
<u>344-08</u>	1 g - 5 kg	-	-	353-480-200	in wooden case
<u>344-084</u>	1 g - 5 kg	-	313-082-400	-	in plastic case
<u>344-09</u>	1 g - 10 kg	-	-	353-490-200	in wooden case
<u>344-026</u>	1 g - 50 g	314-020-600	-	-	in aluminium case
<u>344-036</u>	1 g - 100 g	314-030-600	-	-	in aluminium case
<u>344-046</u>	1 g - 200 g	314-040-600	-	-	in aluminium case
<u>344-056</u>	1 g - 500 g	314-050-600	-	-	in aluminium case
<u>344-066</u>	1 g - 1 kg	314-060-600	-	-	in aluminium case
<u>344-076</u>	1 g - 2 kg	314-070-600	-	-	in aluminium case
<u>344-086</u>	1 g - 5 kg	314-080-600	-	-	in aluminium case
<u>344-096</u>	1 g - 10 kg	314-090-600	-	-	in aluminium case