

KERN 307-01

Individual weight, OIML Class E1, 1 g, Knob, stainless steel polished



Article codes, product groups & Status information

Model article-code	307-01
Model-Serie	307
Assortment	KERN
Product Type	Test weights
Product group	OIML E1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761314184
Predecessor Model	307-1-2.021a
Planned start date	2022-2HJ
Customs tariff number	84239090
Country of Origin	CN

Technical data - Primary

Nominal value	1 g
Tolerance (OIML)	0,01 mg
Case aluminium	317-010-600
Case plastic	317-020-400
Case wood	317-010-100

Technical data - Construction

Material	Stainless steel polished
----------	--------------------------

Technical data - Measuring system

Nominal value	1 g
OIML class	E1
OIML conform	Yes

Functions

Standard

Tolerance (OIML)	0,01 mg
Density (kg/m ³)	8.000
Volume (cm ³)	0,125

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	6×6×6 mm
Gross weight	0,05 kg
Net weight	0,05 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	963-231
-------------------	---------

KERN 307-01

Individual weight, OIML Class E1, 1 g, Knob, stainless steel polished



Option

