

KERN 347-026

Individual weight, OIML Class M1, 2 g, Hook weight, stainless steel fine turned



Article codes, product groups & Status information

Model article-code	347-026
Model-Serie	347-026
Assortment	Scales
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761256163
Predecessor Model	347-426
Planned start date	2022-1HJ
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	2 g
Tolerance (OIML)	1,2 mg
Case plastic	347-030-400

Technical data - Construction

Material	Stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	2 g
OIML class	M1
OIML conform	Yes
Tolerance (OIML)	1,2 mg
Density (kg/m ³)	7.950

Functions

Standard

Volume (cm³) 0,2381

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	6×6×18 mm
Gross weight	0,05 kg
Net weight	0,05 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-632
-------------------	---------

KERN 347-026

Individual weight, OIML Class M1, 2 g, Hook weight, stainless steel fine turned



Option

