

## PRODUCT BROCHURE

---

# KERN 344-08

Date: 09.03.2026

ENGLISH



Weight set, M1, 1 g - 5 kg, Knob, stainless steel fine turned, in wooden case

---

### CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)



Scan here for more  
information

## PRODUCT SHOWCASE



## SPECIFICATIONS

### Article codes, product groups & Status information

|                       |                      |
|-----------------------|----------------------|
| Model article-code    | 344-08               |
| Model-Serie           | 344                  |
| Assortment            | KERN                 |
| Product Type          | Test weights         |
| Product group         | OIML M1 Test weights |
| Product usage type    | Main article         |
| EAN code (Model)      | 4045761008007        |
| Customs tariff number | 84239090             |
| Country of Origin     | DE                   |

### Technical data - Primary

|                 |                |
|-----------------|----------------|
| Nominal value   | 1 g - 5 kg     |
| Case wood       | 353-480-200    |
| Weight set case | in wooden case |

### Technical data - Construction

|                          |                             |
|--------------------------|-----------------------------|
| Material                 | stainless steel fine turned |
| Is weight set            | yes                         |
| Number of weights in set | 17                          |
| Type of case             | in wooden case              |

### Technical data - Measuring system

|               |            |
|---------------|------------|
| Nominal value | 1 g - 5 kg |
| OIML class    | M1         |
| OIML conform  | yes        |

### Technical data - Packing & Shipping

|                              |                |
|------------------------------|----------------|
| Dimensions packaging (W×D×H) | 430×190×200 mm |
| Gross weight                 | 18 kg          |
| Net weight                   | 18 kg          |

|                                               |                |
|-----------------------------------------------|----------------|
| Shipping method                               | Parcel service |
| Packaging component -<br>by weight - plastics | 0 g            |
| Delivery time                                 | 1 d            |

## Optional Services

|                   |         |
|-------------------|---------|
| DAkkS Calibration | 962-621 |
|-------------------|---------|

## FUNCTIONS

|          |
|----------|
| Standard |
|----------|



|        |
|--------|
| Option |
|--------|



## SERVICES

| Model                    | Description                                         |
|--------------------------|-----------------------------------------------------|
| <a href="#">962-621</a>  | DAkkS Calibration KERN 962-621                      |
| <a href="#">962-621R</a> | DAkkS Recalibration KERN 962-621R                   |
| <a href="#">920-716</a>  | Marking (Last digits of serial number) KERN 920-716 |

## SPARE PARTS

| Model                       | Description                     |
|-----------------------------|---------------------------------|
| <a href="#">353-480-200</a> | Box of Weights KERN 353-480-200 |
| <a href="#">315-242</a>     | Tweezers KERN 315-242           |