

PRODUCT BROCHURE

KERN 303-07

Date: 09.03.2026

ENGLISH



Weight set, E1, 1 mg - 2 kg, Knob, stainless steel polished, in wooden case

CONTACT


KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



DESCRIPTION

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	303-07
Model-Serie	303
Assortment	KERN
Product Type	Test weights
Product group	OIML E 1 Test weights
Product usage type	Main article
EAN code (Model)	404576 1002517
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 mg - 2 kg
---------------	-------------

Case wood	313-070-100
Weight set case	in wooden case

Technical data - Construction

Material	stainless steel polished
Is weight set	yes
Number of weights in set	27
Type of case	in wooden case

Technical data - Measuring system

Nominal value	1 mg - 2 kg
OIML class	E1
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	305×190×140 mm
Gross weight	10 kg
Net weight	10 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	963-206
-------------------	---------

FUNCTIONS

Standard



Option



SERVICES

Model	Description
<u>963-206</u>	DAkkS Calibration KERN 963-206
<u>962-206R</u>	DAkkS Recalibration KERN 962-206R
<u>920-715</u>	Marking (Last digits of serial number) KERN 920-715

SPARE PARTS

Model	Description
<u>313-070-100</u>	Box of Weights KERN 313-070-100
<u>315-242</u>	Tweezers KERN 315-242
<u>318-270</u>	Dusting Brush KERN 318-270