

KERN 346-89

Individual weight, OIML Class M1, 50 kg, Block, cast iron lacquered



Article codes, product groups & Status information

Model article-code	346-89
Model-Serie	346
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761241664
Predecessor Model	346-89-2018a
Successor Model	346-89-2024e2; 346-89-2025e
Planned start date	2020-1HJ
Customs tariff number	84239090

Technical data - Primary

Nominal value	50 kg
Tolerance (OIML)	2,5 g

Technical data - Construction

Dimensions	310×155×192 mm
Material	Cast iron lacquered

Technical data - Measuring system

Nominal value	50 kg
OIML class	M1
OIML conform	Yes
Tolerance (OIML)	2,5 g
Density (kg/m ³)	7100

Volume (cm³) 7042,2535

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	330×180×230 mm
Gross weight	50 kg
Net weight	50 kg
Shipping method	Freight
Packaging component - by weight - plastics	0 g
Delivery time	2 d

Technical data - Verification

Verification possible	Yes
-----------------------	-----

Optional Services

DAkkS Calibration	962-646
-------------------	---------

Functions

KERN 346-89

Individual weight, OIML Class M1, 50 kg, Block, cast iron lacquered



Standard

