

KERN OZL 451

Stereo zoom microscope with halogen illumination, for the laboratory, training centres, quality control or agriculture



Article codes, product groups & Status information

Model article-code	OZL 451
Model-Serie	OZL-45
Assortment	Microscopes & Optical Metrology
Product Type	Microscopes
Product group	Stereo microscopes
Product usage type	Main article
EAN code (Model)	4045761170384
Planned start date	2015-1HJ
Customs tariff number	90111000
CE Label	Yes
Country of Origin	CN
Selective distribution	0

Technical data - Primary

Optical System	Greenough
Tube	Binocular
Eyepiece	Eyepiece HSWF 10x / Ø 23 mm with anti-fungus, high eye point
Objectives	0,75x - 5x
Illumination	12 V, 10W Halogen (Durchlicht), 10W Halogen (Auflicht)

Technical data - Construction

System	Upright
Tube angle	45°

Tube rotation	Yes
Stand type	Pillar
Eye distance (max)	75 mm
Eye distance (min)	55 mm
Adjustable dioptre	Both-sided
Focusing mechanics	Coarse drive
Visual field min	5 mm
Visual field max	33 mm
Working distance	113 mm
Sample height max	120 mm
Max resolution	0,004025 mm
Dimensions (W×D×H)	330×270×460 mm

Technical data - Power supply

Power supply kind	Built-in power supply
Input Voltage	220 V
Power Plug	Cold device cable
Protection Class	Protection class I - protective conductor
Cable attachment	Removable

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Storage Temp Max	40 °C

Technical data - Packing & Shipping

KERN OZL 451

Stereo zoom microscope with halogen illumination, for the laboratory, training centres, quality control or agriculture



Dimensions packaging (W×D×H)	340×340×480 mm
Gross weight	6 kg
Net weight	4,6 kg
Shipping method	Parcel service
Packaging component - by weight - paper	0 g
Packaging component - by weight - cardboard	749,95 g
Packaging component - by weight - plastics	9,905 g
Packaging component - by weight - styrofoam	655,145 g
Packaging component - by weight - foam	0 g
Delivery time	1 d

Technical data - Verification

WEEE	Yes
------	-----

Technical data - Illumination

Illumination directions	Reflected light; Transmitted light
Transmitted light type	Halogen
Transmitted light intensity	10 W
Incident light type	Halogen
Incident light intensity	10 W
Dimmable	Incident light + Transmitted light, collectively

Functions

Standard

