

PRODUCT BROCHURE

KERN 357-08

Date: 10.03.2026

ENGLISH



Individual weight, OIML Class M2, 200 g, Knob, stainless steel fine turned

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	357-08
Model-Serie	357-08
Assortment	KERN
Product Type	Test weights
Product group	OIML M2 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761289000
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	200 g
Tolerance (OIML)	30 mg
Case aluminium	317-080-600
Case plastic	347-080-400
Case wood	337-080-200

Technical data - Construction

Material	stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	200 g
OIML class	M2

OIML conform	yes
Tolerance (OIML)	30 mg
Density (kg/m ³)	7950
Volume (cm ³)	23,8095

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	31×31×56 mm
Gross weight	0.20 kg
Net weight	0.20 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-638
-------------------	---------

FUNCTIONS

Standard



Option



ACCESSORIES

Model	Description
<u>337-080-200</u>	Wooden Box for Set of Weights KERN 337-080-200
<u>347-080-400</u>	Plastic Weight Box KERN 347-080-400
<u>317-080-600</u>	Aluminium Protective Box KERN 317-080-600
<u>315-242</u>	Tweezers KERN 315-242

SERVICES

Model	Description
<u>962-638</u>	DAkkS Calibration KERN 962-638
<u>962-638R</u>	DAkkS Recalibration KERN 962-638R