

## PRODUCT BROCHURE

---

# KERN SET-52

Date: 28.04.2025

ENGLISH

# Software EasyTouch SET-52 Inventory

EasyTouch Inventory – stock-taking: Efficient and correct

---

### CONTACT





KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)



Scan here for more  
information

## PRODUCT SHOWCASE



## DESCRIPTION

- Note: Prerequisite for this set are the basic program SET-01 Base and SET-32 Count
- ET Inventory: SET-52 Inventory is an assistance program and database, which enables users without any previous knowledge to carry out stock-take processes quickly and without errors. Users concentrate solely on identifying the stock-take objects, e.g. by scanning the barcode, while the program handles the entire counting process independently
- Simple process: The process is remarkably easy: Scan the barcode, place the object on the balance and you're done. Then go straight on to the next item
- Organising the stock-take: Each stock-take can be divided into different "jobs" – depending on the storage location, product group, etc. Each stock-taking job is carried out, completed and saved individually. This makes it possible to divide large stock-takes into meaningful subsections and have them processed by different teams
- Authorisation process: Stock-takes can include an authorisation process so that the results have to be authorised by an authorised user role when the process has been completed
- Stock list: EasyTouch outputs a clear list of all recorded items with the counted stock for each item. It is possible, for example, to compare a current stock-take with the stock-take from the previous year. The system

then calculates the stock differences between the two stock-takes automatically

- ERP Connect: With this optional extension, stock lists can be imported from existing order management systems for example and results, including current stock levels, can be exported back to the system
- Weighing and counting: The stock value can either be taken from the balance as a count result, can be entered by hand or taken from a comparison stock-take
- Data import: The master data for the stock-take item can be imported into EasyTouch either from an ERP system using an interface, or from an Excel list

## SPECIFICATIONS

### Article codes, product groups & Status information

Model article-code	SET-52
Model-Serie	ETS
Assortment	KERN
Product Type	Software
Product group	EasyTouch
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761486362
Customs tariff number	85234910

### Technical data - Packing & Shipping

Delivery time	1 d
---------------	-----

## FUNCTIONS

### Standard



### Option



# SOFTWARES

Model	Description
SET-14	KERN EasyTouch SET-14 Individual Print