

KERN FOB 7K-4NLO

Compact organ scale made of stainless steel, particularly well protected and hygienic to clean thanks to IP67 protection and smooth stainless steel surfaces



Article codes, product groups & Status information

Model article-code	FOB 7K-4NLO
Model-Serie	FOB
Assortment	KERN
Product Type	Medical products
Product group	Personal scales
Product usage type	Main article
EAN code (Model)	4045761264366
Planned start date	2021-2HJ
Customs tariff number	84238129
CE Label	Yes
Country of Origin	CN

Technical data - Primary

Weighing capacity	5.000 g; 7.500 g
Readability	0,5 g; 1 g

Technical data - Construction

Material weighing plate	Stainless steel
Weighing surface (d)	252×200 mm
Dimensions of weighing plate (W×D)	252×200 mm
Material casing	Stainless steel
IP protection indicator	IP67
IP protection platform	IP67

Technical data - Display

Display height	25 mm
----------------	-------

Technical data - Power supply

Battery	4×15 V AA
Operating time	48 h
Input Voltage	230 V AC, 50 Hz

Technical data - Measuring system

Calibration method	Adjusting with external weight
Linearity	± 1,5 g; 3 g
Recommended adjusting weight	5 kg (M1)
Repeatability	0,5 g; 1 g
Weighing units	G;lb
Number of weighing ranges	2
Weighing system	Strain gauge

Technical data - Environmental conditions

Maximum humidity	95 %
Maximum operating temperature	35 °C
Minimum ambient temperature	5 °C
Minimum humidity	25 %

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	415×305×165 mm
Gross weight	4,6 kg
Net weight	3,8 kg
Shipping method	Parcel service

KERN FOB 7K-4NLO

Compact organ scale made of stainless steel, particularly well protected and hygienic to clean thanks to IP67 protection and smooth stainless steel surfaces



Packaging component - by weight - paper 100 g

Packaging component - by weight - cardboard 347,5 g

Packaging component - by weight - plastics 347,5 g

Packaging component - by weight - styrofoam 0 g

Packaging component - by weight - glass 0 g

Packaging component - by weight - wood 0 g

Packaging component - by weight - foam 0 g

Packaging component - by weight - metal 0 g

Delivery time 1 d

Optional Services

DAkkS Calibration 963-128

Functions

Standard



Option

