

PRODUCT BROCHURE

KERN 343-08

Date: 09.03.2026

ENGLISH



Weight set, M1, 1 mg - 5 kg, Knob, stainless steel fine turned, in wooden case


CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information



ENGLISH

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

Article codes, product groups & Status information

Model article-code	343-08
Model-Serie	343
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main article
EAN code (Model)	4045761007581
Customs tariff number	84239090
Country of Origin	DE

Nominal value	1 mg - 5 kg
---------------	-------------

Case wood	353-480-200
Weight set case	in wooden case

Material	stainless steel fine turned
Is weight set	yes
Number of weights in set	28
Type of case	in wooden case

Nominal value	1 mg - 5 kg
OIML class	M1
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	430×190×200 mm
Gross weight	18 kg
Net weight	18 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-607
-------------------	---------

FUNCTIONS

Standard



Option



SERVICES

Model	Description
<u>962-607</u>	DAkkS Calibration KERN 962-607
<u>962-607R</u>	DAkkS Recalibration KERN 962-607R
<u>920-708</u>	Marking (Last digits of serial number) KERN 920-708

SPARE PARTS

Model	Description
<u>353-480-200</u>	Box of Weights KERN 353-480-200
<u>315-242</u>	Tweezers KERN 315-242
<u>318-270</u>	Dusting Brush KERN 318-270