


KERN 346-83

Individual weight, OIML Class M1, 500 kg, Heavy Duty Weight, steel, lacquered



Article codes, product groups & Status information

Model article-code	346-83
Model-Serie	346
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761479579
Predecessor Model	346-83-2024a
Planned start date	2025-2HJ
Customs tariff number	84239090

Technical data - Primary

Nominal value	500 kg
Tolerance (OIML)	25 g

Technical data - Construction

Dimensions (W×D×H)	760×500×320 mm
Material	Steel, lacquered

Technical data - Measuring system

Nominal value	500 kg
OIML class	M1
OIML conform	Yes
Tolerance (OIML)	25 g
Density (kg/m ³)	7100
Volume (cm ³)	69444

Functions

Standard

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	760×500×320 mm
Gross weight	500 kg
Net weight	500 kg
Shipping method	Freight
Packaging component - by weight - plastics	0 g
Delivery time	2 d

Optional Services

DAkkS Calibration	962-693
-------------------	---------

KERN 346-83

Individual weight, OIML Class M1, 500 kg, Heavy Duty Weight, steel, lacquered

