

PRODUCT BROCHURE

KERN 325-096

Date: 10.03.2026

ENGLISH



Weight set, F1, 1 mg - 10 kg, Eco Shape, stainless steel polished, in aluminium case

CONTACT


KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



ENGLISH

DESCRIPTION

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	325-096
Model-Serie	325
Assortment	KERN
Product Type	Test weights
Product group	OIML F1 Test weights
Product usage type	Main article
EAN code (Model)	404576 1224 155
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	1 mg - 10 kg
---------------	--------------

Case aluminium	313-090-600
Weight set case	in aluminium case

Technical data - Construction

Material	stainless steel polished
Is weight set	yes

Technical data - Measuring system

OIML class	F1
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	610×260×235 mm
Gross weight	28 kg

Net weight	20 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	8 kg
Delivery time	1 d

Optional Services

DAkkS Calibration	962-408
-------------------	---------

FUNCTIONS

Standard



Option



SERVICES

Model	Description
<u>962-408</u>	DAkkS Calibration KERN 962-408
<u>962-408R</u>	DAkkS Recalibration KERN 962-408R
<u>921-408</u>	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-408
<u>920-709</u>	Marking (Last digits of serial number) KERN 920-709

SPARE PARTS

Model	Description
<u>313-090-600</u>	Box of Weights KERN 313-090-600