

KERN 347-077

Individual weight, OIML Class M1, 100 g, Slot weight (Type B), stainless steel fine turned



Article codes, product groups & Status information

Model article-code	347-077
Model-Serie	347-077
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761366695
Planned start date	2024-1HJ
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	100 g
Tolerance (OIML)	5 mg

Technical data - Construction

Material	Stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	100 g
OIML class	M1
OIML conform	Yes
Tolerance (OIML)	5 mg
Density (kg/m ³)	8.000
Volume (cm ³)	11,9048

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	60×60×8,1 mm
Gross weight	0,10 kg
Net weight	0,10 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-637
-------------------	---------

Functions

KERN 347-077

Individual weight, OIML Class M1, 100 g, Slot weight (Type B), stainless steel fine turned



Standard



Option

