

SAUTER HO-A23

Date: 29.01.2025

ENGLISH



Hardness test block for stainless steel HO; (58~63)HRC without certificate

CONTACT

 Address

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

 Phone and fax

Phone : +49 - [0] 7433 - 9933 - 0
Fax : +49 - [0] 7433 - 9933 - 149

 Internet

Email : info@kern-sohn.com
Web : www.kern-sohn.com

DESCRIPTION

Hardness comparison plate HO-A23

SPECIFICATION

Article codes, product groups & Status information

Model article-code	HO-A23
Model-Serie	HO
Assortment	Sauter
Product Type	Hardness measurement
Product group	Hardness testers (UCI)
Product usage type	Spare part
EAN code (Model)	404576 125 1625
Customs tariff number	90248000

Technical data - Construction

Dimensions of housing	90×90×16 mm
Material housing	stainless steel

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	150×150×20 mm
Gross weight	0.90 kg
Net weight	0.85 kg
Shipping method	Parcel service
UN number	No dangerous goods
Packaging component - by weight - plastics	10 g
Delivery time	1 d

RELATED ARTICLES DATA

Model name	Measuring range Rockwell - scale A	Measuring range Rockwell - scale B	Measuring range Rockwell - scale C
HO-A22			
RHO-7001			
HO-A23			
HO-A03			