

PRODUCT BROCHURE

KERN SCD-4.0-PRO-DL

Date: 18.03.2025

ENGLISH



Flexible recording or transfer of measurements, particularly to Microsoft® Excel or Access

CONTACT

KERN & SOHN GmbH

Ziegelei 1

72336 Balingen

Germany

Phone : +49 7433 9933-0

Fax : +49 7433 9933-149

Email : info@kern-sohn.comWeb : www.kern-sohn.com

- Many different transfer and recording options:
 - Microsoft Excel/Access/Word in the foreground or background
 - Other Windows applications (through key simulation), e.g. shipping software or ERP system (SAP, Sage, etc.)
 - File recording (e.g. as text or CSV file)
 - Print out on text or label printer
- Additional features of the BalanceConnection PRO version::
 - ODBC/SQL databases such as SQL Server or MySQL
 - Transfer to HTTP Webservices/Web forms
 - Histogram
 - Enables the connection of medical balances to practice EDP systems using the device data transfer protocol (GDT) and also the HL7 protocol
 - Implementation of programs/scripts using freely-definable command lines
 - Command sequences
 - Defining the conditions for specific events and the reaction to these events
 - Filter for stability recognition (configurable value storage)
 - Central configuration repository, e.g. on the network
 - On-screen displays with configurable user interface for flexible balance operation
- Benefits
 - Ideal software where weighing data which has been collected needs to be transferred to a Windows PC
 - Particularly suitable for professional integration of weighing systems into existing networks where there is a high demand for automated data transfer and further processing, for integration of databases
 - Automatic recording, e.g. for long-term tests
 - GLP-compliant lab software
 - For each license, you can use as many balances as you need at one PC/workstation
 - Can be used across different branches and applications
- Scope of supplies: Download link for 1 license

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	SCD-4.0-PRO-DL
Model-Serie	SCD
Assortment	KERN
Product Type	Software
Product group	BalanceConnection

Product usage type	Main item;Accessories/Services
EAN code (Model)	404576 1222489
Predecessor Model	SCD-4.0; SCD-4.0-PRO; SCD-4.0-PROS05; SCD-4.0S05
Customs tariff number	852349 10

Technical data - Packing & Shipping

Delivery time 1 d

FUNCTIONS

Standard

