

# KERN 347-047

Individual weight, OIML Class M1, 10 g, Slot weight (Type B), stainless steel fine turned



## Article codes, product groups & Status information

|                       |                                |
|-----------------------|--------------------------------|
| Model article-code    | 347-047                        |
| Model-Serie           | 347-047                        |
| Assortment            | KERN                           |
| Product Type          | Test weights                   |
| Product group         | OIML M1 Test weights           |
| Product usage type    | Main item;Accessories/Services |
| EAN code (Model)      | 4045761366633                  |
| Planned start date    | 2024-1HJ                       |
| Customs tariff number | 84239090                       |
| Country of Origin     | DE                             |

## Technical data - Primary

|                  |      |
|------------------|------|
| Nominal value    | 10 g |
| Tolerance (OIML) | 2 mg |

## Technical data - Construction

|          |                             |
|----------|-----------------------------|
| Material | Stainless steel fine turned |
|----------|-----------------------------|

## Technical data - Measuring system

|                              |        |
|------------------------------|--------|
| Nominal value                | 10 g   |
| OIML class                   | M1     |
| OIML conform                 | Yes    |
| Tolerance (OIML)             | 2 mg   |
| Density (kg/m <sup>3</sup> ) | 8.000  |
| Volume (cm <sup>3</sup> )    | 1,1905 |

## Technical data - Packing & Shipping

|                              |            |
|------------------------------|------------|
| Dimensions packaging (W×D×H) | 28×28×6 mm |
|------------------------------|------------|

|              |         |
|--------------|---------|
| Gross weight | 0,05 kg |
|--------------|---------|

|            |         |
|------------|---------|
| Net weight | 0,05 kg |
|------------|---------|

|                 |                |
|-----------------|----------------|
| Shipping method | Parcel service |
|-----------------|----------------|

|                                            |     |
|--------------------------------------------|-----|
| Packaging component - by weight - plastics | 0 g |
|--------------------------------------------|-----|

|               |     |
|---------------|-----|
| Delivery time | 1 d |
|---------------|-----|

## Optional Services

|                   |         |
|-------------------|---------|
| DAkkS Calibration | 962-634 |
|-------------------|---------|

## Functions

# KERN 347-047

Individual weight, OIML Class M1, 10 g, Slot weight (Type B), stainless steel fine turned



## Standard



## Option

