

PRODUCT BROCHURE

KERN 347-095-100

Date: 24.03.2025

ENGLISH



Supporting bar, OIML Class M1, 500 g, stainless steel fine turned

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com

DESCRIPTION

- for fixing slotted weights
- Own weight beam bar 500 g
- Maximum total load 20 kg (is exclusive of the own weight of the beam bar, e.g. the maximum possible total weight is calculated from
“Maximum total load” + “own weight beam bar”)
- Largest slotted weight possible 10 kg

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	347-095-100
Model-Serie	347-095-100
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	404576 1256354
Predecessor Model	347-495-100
Customs tariff number	84239090

Technical data - Primary

Nominal value	500 g
Tolerance (OIML)	25 mg
Case plastic	347-515-400

Technical data - Construction

Material	stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	500 g
OIML class	M1
OIML conform	yes
Tolerance (OIML)	25 mg
Density (kg/m ³)	8000
Volume (cm ³)	59,5238

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	40×40×663 mm
Gross weight	0.50 kg
Net weight	0.50 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-639
-------------------	---------

FUNCTIONS

Standard



ACCESSORIES

Model	Description
347-515-400	Plastic Weight Box KERN 347-515-400
347-015	Test Weight KERN 347-015
347-025	Test Weight KERN 347-025
347-035	Test Weight KERN 347-035
347-045	Test Weight KERN 347-045
347-055	Test Weight KERN 347-055
347-065	Test Weight KERN 347-065
347-075	Test Weight KERN 347-075
347-085	Test Weight KERN 347-085
347-095	Test Weight KERN 347-095
347-115	Test Weight KERN 347-115
347-125	Test Weight KERN 347-125
347-135	Test Weight KERN 347-135
347-145	Test Weight KERN 347-145

SERVICES

Model	Description
962-639	DAkkS Calibration KERN 962-639
962-639R	DAkkS Recalibration KERN 962-639R
920-32	Marking by etching KERN 920-32

RELATED ARTICLES DATA

Model name	Nominal value	Tolerance (OIML)	Case plastic
347-095-100	500 g	25 mg	347-515-400