


KERN 346-022-005

Weight container for rectangular weights, calibrated to OIML class M1



Article codes, product groups & Status information

Model article-code	346-022-005
Model-Serie	346
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main article
EAN code (Model)	4045761225053
Planned start date	2019-1HJ
Customs tariff number	84239090

Net weight	20 kg
Shipping method	Freight
Packaging component - by weight - plastics	35 kg
Delivery time	2 d

Optional Services

DAkkS Calibration	962-645
-------------------	---------

Technical data - Primary

Nominal value	20 kg
Tolerance (OIML)	1 g

Technical data - Construction

Dimensions (W×D×H)	730×300×330 mm
Material	Stainless steel

Technical data - Measuring system

Nominal value	20 kg
OIML class	M1
OIML conform	Yes
Tolerance (OIML)	1 mg

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	730×300×330 mm
Gross weight	55 kg

Functions

KERN 346-022-005

Weight container for rectangular weights, calibrated to OIML class M1



Standard

