

KERN 347-015-N00

Individual weight (Newton), 1 cN, Slot weight (Type A), stainless steel fine turned



Article codes, product groups & Status information

Model article-code	347-015-N00
Model-Serie	347-015-N00
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761277137
Predecessor Model	347-415-N00
Planned start date	2022-1HJ
Customs tariff number	84239090

Technical data - Primary

Tolerance (OIML)	1 mg
Case plastic	347-030-400

Technical data - Construction

Material	Stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

OIML class	M1
OIML conform	No
Tolerance (OIML)	1 mg
Density (kg/m ³)	7950
Volume (cm ³)	0,119

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	13×13×2 mm
------------------------------	------------

Gross weight	0.05 kg
Net weight	0.05 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-631N
-------------------	----------

Functions

KERN 347-015-N00

Individual weight (Newton), 1 cN, Slot weight (Type A), stainless steel fine turned



Standard

