

# KERN 327-11

Individual weight, OIML Class F1, 1 kg, Knob, stainless steel polished



## Article codes, product groups & Status information

|                       |                                |
|-----------------------|--------------------------------|
| Model article-code    | 327-11                         |
| Model-Serie           | 327                            |
| Assortment            | KERN                           |
| Product Type          | Test weights                   |
| Product group         | OIML F1 Test weights           |
| Product usage type    | Main item;Accessories/Services |
| EAN code (Model)      | 4045761004818                  |
| Successor Model       | 327-11-2025e                   |
| Planned start date    | 2004-1HJ                       |
| Customs tariff number | 84239090                       |

## Technical data - Primary

|                  |             |
|------------------|-------------|
| Nominal value    | 1 kg        |
| Tolerance (OIML) | 5 mg        |
| Case aluminium   | 317-110-600 |
| Case plastic     | 347-110-400 |
| Case wood        | 317-110-100 |

## Technical data - Construction

|          |                          |
|----------|--------------------------|
| Material | Stainless steel polished |
|----------|--------------------------|

## Technical data - Measuring system

|                  |      |
|------------------|------|
| Nominal value    | 1 kg |
| OIML class       | F1   |
| OIML conform     | Yes  |
| Tolerance (OIML) | 5 mg |

## Functions

Standard

Density (kg/m<sup>3</sup>) 7950

Volume (cm<sup>3</sup>) 125

## Technical data - Packing & Shipping

|                                            |                |
|--------------------------------------------|----------------|
| Dimensions packaging (W×D×H)               | 48×48×83 mm    |
| Gross weight                               | 1.00 kg        |
| Net weight                                 | 1.00 kg        |
| Shipping method                            | Parcel service |
| Packaging component - by weight - plastics | 0 g            |
| Delivery time                              | 1 d            |

## Technical data - Verification

Verification possible Yes

## Optional Services

DAkkS Calibration 962-441

# KERN 327-11

Individual weight, OIML Class F1, 1 kg, Knob, stainless steel polished

