

PRODUCT BROCHURE

KERN 313-06

Date: 09.03.2026

ENGLISH



Weight set, E2, 1 mg - 1 kg, Knob, stainless steel polished, in wooden case


CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information



ENGLISH

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

Article codes, product groups & Status information

Case wood	313-060-100
Weight set case	in wooden case

Material	stainless steel polished
Is weight set	yes
Number of weights in set	25
Type of case	in wooden case

Nominal value	1 mg - 1 kg
OIML class	E2
OIML conform	yes

Nominal value	1 mg - 1 kg
---------------	-------------

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	230×180×110 mm
Gross weight	4.4 kg
Net weight	4.4 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-305
-------------------	---------

FUNCTIONS

Standard



Option



SERVICES

Model	Description
<u>962-305</u>	DAkkS Calibration KERN 962-305
<u>962-305R</u>	DAkkS Recalibration KERN 962-305R
<u>921-305</u>	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-305
<u>920-706</u>	Marking (Last digits of serial number) KERN 920-706

SPARE PARTS

Model	Description
<u>313-060-100</u>	Box of Weights KERN 313-060-100
<u>315-242</u>	Tweezers KERN 315-242
<u>318-270</u>	Dusting Brush KERN 318-270