

PRODUCT BROCHURE

OPTICS OZB-A5 130

Date: 30.04.2025

ENGLISH



Individuality, variety and flexible working through our modular construction system ? ECO-Universal stands

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

Stereomicroscope Stand OZB-A5130

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	OZB-A5130
Model-Serie	OZB-S
Assortment	Optics
Product Type	Microscopes
Product group	Stereo microscopes
Product usage type	Spare part;Accessories/Services
EAN code (Model)	4045761232136
Predecessor Model	OZB-A5103
Customs tariff number	90119000
CE Label	yes

Technical data - Construction

Stand type	pillar form
Stand area	283×292 mm
Holder diameter	76 mm
Focusing mechanics	coarse drive
Dimensions (W×D×H)	283×292×303,5 mm

Technical data - Power supply

Protection Class	Class III - Low voltage
------------------	-------------------------

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Stotage Temp Max	40 °C

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	385×370×330 mm
Gross weight	4.2 kg
Net weight	3 kg
Shipping method	Parcel service
Packaging component - by weight - paper	0 g
Packaging component - by weight - cardboard	780.168 g
Packaging component - by weight - plastics	14.8368 g
Packaging component - by weight - styrofoam	441.395 g
Packaging component - by weight - foam	0 g
Delivery time	1 d

Technical data - Verification

WEEE	yes
------	-----

Technical data - Illumination

Illumination directions	Reflected light;Transmitted light
-------------------------	-----------------------------------

Transmitted light type LED

Transmitted light
intensity 3W

Incident light type LED

Incident light intensity 3W

FUNCTIONS

Standard



SPARE PARTS

Model	Description
OZB-A5807	Power Adapter KERN OZB-A5807
OZB-A5808	Power Adapter KERN OZB-A5808