

KERN 326-076

Weight set, F1, 1 g - 2 kg, Eco Shape, stainless steel polished, in aluminium case



Article codes, product groups & Status information

Model article-code	326-076
Model-Serie	326
Assortment	KERN
Product Type	Test weights
Product group	OIML F1 Test weights
Product usage type	Main article
EAN code (Model)	4045761224094
Planned start date	2019-1HJ
Customs tariff number	84239090

Technical data - Primary

Nominal value	1 g - 2 kg
Case aluminium	314-070-600
Weight set case	In aluminium case

Technical data - Construction

Material	Stainless steel polished
Is weight set	Yes
Type of case	In aluminium case

Technical data - Measuring system

Nominal value	1 g - 2 kg
OIML class	F1
OIML conform	Yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	278×180×138 mm
------------------------------	----------------

Functions

Gross weight	6 kg
Net weight	1.4 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	4290 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-420
-------------------	---------

KERN 326-076

Weight set, F1, 1 g - 2 kg, Eco Shape, stainless steel polished, in aluminium case



Standard

