

PRODUCT BROCHURE

KERN 347-136

Date: 10.03.2026

ENGLISH



Individual weight, OIML Class M1, 5 kg, Hook weight, stainless steel fine turned

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

- Individual weight, OIML Class M1, 5 kg, Hook weight, stainless steelfine turned (OIML)

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	347-136
Model-Serie	347-136
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761255418
Predecessor Model	347-536
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	5 kg
Tolerance (OIML)	250 mg
Case plastic	347-140-400

Technical data - Construction

Material	stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	5 kg
OIML class	M1
OIML conform	yes
Tolerance (OIML)	250 mg
Density (kg/m ³)	8000
Volume (cm ³)	595,2381

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	80×80×130 mm
Gross weight	5 kg
Net weight	5 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-643
-------------------	---------

FUNCTIONS

Standard



Option



ACCESSORIES

Model	Description
<u>317-280</u>	Gloves KERN 317-280
<u>317-290</u>	Gloves KERN 317-290
<u>347-140-400</u>	Plastic Weight Box KERN 347-140-400

SERVICES

Model	Description
<u>962-643</u>	DAkkS Calibration KERN 962-643
<u>962-643R</u>	DAkkS Recalibration KERN 962-643R