

PRODUCT BROCHURE

KERN YPS-03

Date: 16.03.2025

ENGLISH



Antivibration table to absorb vibrations and oscillations, which would otherwise distort the weighing result

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com

DESCRIPTION

- The KERN YPS-03 antivibration table has been constructed to absorb vibrations and oscillations, which would otherwise distort the weighing result
- The YPS-03 antivibration table consists of:
 - a heavy, polished granite slab (24 kg) and is positioned on absorbing rubber components and is used as a surface on which to stand the analytical balance
 - a work table, the work surface of which surrounds the slab and offers many storage options. The sturdy powder-coated steel frame has four height-adjustable feet, which can be used to accurately level the antivibration table if the floor is at all uneven
- You can work in an ergonomic way at the antivibration table. The inlet on the front of the table offers enough space for you to sit down and operate the scale without touching the antivibration table
- All surfaces are hygienic and easy to clean
- Suitable for all KERN analytical balances with overall dimensions =W×D 400×450 mm
- Overall dimensions W×D×H 800×600×785 mm (height is levellable, approx. 20 mm)
- Dimensions of granite slab, W×D×H 410×270×70 mm
- Net weight approx. 46 kg

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	YPS-03
Model-Serie	YPS
Assortment	KERN
Product Type	General weighing accessories
Product group	Accessories (printer, etc.)
Product usage type	Accessories/DL
EAN code (Model)	404576 1246898
Predecessor Model	YPS-03-2019a
Customs tariff number	84239010

Technical data - Construction

Overall dimensions mounted (W×D×H)	800×600×785 mm
------------------------------------	----------------

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	800×600×300 mm
Gross weight	55 kg
Net weight	46 kg
Shipping method	Freight
Packaging component - by weight - cardboard	3750 g
Packaging component - by weight - plastics	3750 g
Delivery time	2 d

FUNCTIONS

Standard



SPARE PARTS

Model	Description
RTYPS-9001	Adjustable base foot KERN RTYPS-9001
RTYPS-8002	Rubber shock absorber KERN RTYPS-8002