

KERN ODC 895

The digital USB microscope for rapid testing or for hobby use



Article codes, product groups & Status information

Model article-code	ODC 895
Model-Serie	ODC-89
Assortment	Optics
Product Type	Microscopes
Product group	Microscope cameras
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761190566
Predecessor Model	ODC 893; ODC 894
Planned start date	2016-1HJ
Customs tariff number	90118000
CE Label	Yes
Country of Origin	CN

Technical data - Primary

Camera type	USB Digital Microscope
Camera resolution	2 MP
Sensor type	CMOS
Sensor size	1/3,2"

Technical data - Construction

Color capability	Colour
Interfaces	USB 2.0
Operating system	Win 10;Win 11;Win 8
Applicable software	Microscope VIS 2.0 Pro (OXM 902);Windows Kamera App

Technical data - Verification

WEEE	Yes
------	-----

Dimensions (W×D×H) 80×150×180 mm

Technical data - Power supply

Power Plug	USB A
Protection Class	Class III - Low voltage
Cable attachment	Permanently connected

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	160×158×85 mm
Gross weight	1,00 kg
Net weight	0,90 kg
Shipping method	Parcel service
Packaging component - by weight - paper	0 g
Packaging component - by weight - cardboard	131,984 g
Packaging component - by weight - plastics	438 mg
Packaging component - by weight - styrofoam	13,432 g
Packaging component - by weight - foam	0 g
Delivery time	1 d

Technical data - Environmental conditions

Storage Temp Min	-5 °C
Storage Temp Max	40 °C

KERN ODC 895

The digital USB microscope for rapid testing or for hobby use



Functions

Standard

