

PRODUCT BROCHURE

KERN ABP-A01

Date: 18.03.2025

ENGLISH



Ioniser to neutralise electrostatic charge

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com

PRODUCT SHOWCASE



DESCRIPTION

- Goods to be weighed which are non-conductive, such as plastic, china, glass, etc often carry electrostatic charge. The electromagnetic field which then occurs between the goods to be weighed and the balance can cause the weighing result to be distorted by up to a gramme. Therefore, when weighing in milligrammes or even smaller increments, we recommend neutralising the electrostatic charge with the ioniser KERN ABP-A01
- It can be conveniently placed next to the scale by means of the stand or directed by hand from several sides to a critical sample to be weighed

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	ABP-A01
Model-Serie	ABP
Assortment	KERN
Product Type	Laboratory balances
Product group	Accessories (printer, etc.)
Product usage type	Main item;Accessories/Services
EAN code (Model)	404576 12 13975
Customs tariff number	90160090

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	325×263×158 mm
Gross weight	1.00 kg
Net weight	0.85 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	75 g
Packaging component - by weight - plastics	75 g
Packaging component - by weight - glass	0 g
Packaging component - by weight - wood	0 g

Delivery time 1 d

FUNCTIONS

Standard



SERVICES

Model	Description
970-109	Warranty extension KERN (+2 years) 970-109