

PRODUCT BROCHURE

KERN 347-137-100-N00

Date: 23.04.2026

ENGLISH



Supporting bar (Newton), 50 N, stainless steel fine turned

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

- for fixing slotted weights
- A precise curvature at the lower end of the beam bar ensures that the slotted weights rest securely. This design prevents the weights from slipping and ensures stability – even in slight inclined positions

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	347-137-100-N00
Model-Serie	347-7-N00
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761366824
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Tolerance (OIML)	250 mg
------------------	--------

Technical data - Construction

Material	stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

OIML class	M1
Density (kg/m ³)	8000

Volume (cm ³)	595,2381
---------------------------	----------

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	138×138×691.5 mm
Gross weight	5 kg
Net weight	5 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	6 d

Optional Services

DAkkS Calibration	962-643N
-------------------	----------

FUNCTIONS

Standard



Option



ACCESSORIES

Model	Description
<u>317-280</u>	Gloves KERN 317-280
<u>317-290</u>	Gloves KERN 317-290
<u>347-037-N00</u>	Test Weight KERN 347-037-N00
<u>347-047-N00</u>	Test Weight KERN 347-047-N00
<u>347-057-N00</u>	Test Weight KERN 347-057-N00
<u>347-067-N00</u>	Test Weight KERN 347-067-N00
<u>347-147-N00</u>	Test Weight KERN 347-147-N00
<u>347-137-N00</u>	Test Weight KERN 347-137-N00
<u>347-127-N00</u>	Test Weight KERN 347-127-N00
<u>347-117-N00</u>	Test Weight KERN 347-117-N00
<u>347-097-N00</u>	Test Weight KERN 347-097-N00
<u>347-087-N00</u>	Test Weight KERN 347-087-N00
<u>347-077-N00</u>	Test Weight KERN 347-077-N00

SERVICES

Model	Description
<u>962-643RN</u>	DAkkS Recalibration KERN 962-643RN
<u>962-643N</u>	DAkkS Calibration KERN 962-643N

RELATED ARTICLES DATA

Model name	Tolerance (OIML)
<u>347-037-N00</u>	1,6 mg
<u>347-047-N00</u>	2 mg
<u>347-057-N00</u>	2,5 mg
<u>347-067-N00</u>	3 mg
<u>347-147-N00</u>	500 mg
<u>347-147-100-N00</u>	500 mg
<u>347-137-N00</u>	250 mg
<u>347-137-100-N00</u>	250 mg

Model name	Tolerance (OIML)
<u>347-127-N00</u>	100 mg
<u>347-117-N00</u>	50 mg
<u>347-117-100-N00</u>	50 mg
<u>347-097-N00</u>	25 mg
<u>347-087-N00</u>	10 mg
<u>347-077-N00</u>	5 mg