

PRODUCT BROCHURE

KERN 347-087

Date: 10.03.2026

ENGLISH



Individual weight, OIML Class M1, 200 g, Slot weight (Type B), stainless steel fine turned

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

- Individual weight, OIML Class M1, 200 g, Slot weight (Type B), stainless steelfine turned (OIML)

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	347-087
Model-Serie	347-087
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761366718
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	200 g
Tolerance (OIML)	10 mg

Technical data - Construction

Material	stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	200 g
OIML class	M1

OIML conform	yes
Tolerance (OIML)	10 mg
Density (kg/m ³)	8000
Volume (cm ³)	23,8095

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	60×60×14.1 mm
Gross weight	0.20 kg
Net weight	0.20 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	962-638
-------------------	---------

FUNCTIONS

Standard



Option



ACCESSORIES

Model	Description
<u>347-137-100</u>	Supporting Bar KERN 347-137-100
<u>347-147-100</u>	Supporting Bar KERN 347-147-100
<u>347-117-100</u>	Supporting Bar KERN 347-117-100
<u>315-242</u>	Tweezers KERN 315-242

SERVICES

Model	Description
<u>962-638</u>	DAkkS Calibration KERN 962-638
<u>962-638R</u>	DAkkS Recalibration KERN 962-638R