

PRODUCT BROCHURE

KERN 307-16

Date: 09.03.2026

ENGLISH



Individual weight, OIML Class E1, 50 kg, Knob, stainless steel polished

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

- Individual weight, OIML Class E1, 50 kg, Knob, Stainless steel polished
- Note: For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	307-16
Model-Serie	307
Assortment	KERN
Product Type	Test weights
Product group	OIML E1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761002807
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary

Nominal value	50 kg
Tolerance (OIML)	25 mg
Case aluminium	317-160-600
Case wood	317-160-100

Technical data - Construction

Material	stainless steel polished
----------	--------------------------

Technical data - Measuring system

Nominal value	50 kg
OIML class	E1
OIML conform	yes
Tolerance (OIML)	25 mg
Density (kg/m ³)	8000
Volume (cm ³)	6250

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	178×178×298 mm
Gross weight	50 kg
Net weight	50 kg
Shipping method	Freight
Packaging component - by weight - plastics	0 g
Delivery time	30 d

Technical data - Verification

Verification possible	yes
-----------------------	-----

Optional Services

DAkkS Calibration	963-246
-------------------	---------

FUNCTIONS

Standard



Option



ACCESSORIES

Model	Description
<u>317-290</u>	Gloves KERN 317-290
<u>317-160-100</u>	Wooden Box for Set of Weights KERN 317-160-100
<u>317-160-600</u>	Aluminium Protective Box KERN 317-160-600

SERVICES

Model	Description
<u>963-246</u>	DAkkS Calibration KERN 963-246
<u>962-246R</u>	DAkkS Recalibration KERN 962-246R
<u>921-246</u>	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-246