

KERN 315-211-407

Safety Set 100 g (M1); 1000 g (F1)



Article codes, product groups & Status information

Model article-code	315-211-407
Model-Serie	315
Assortment	KERN
Product Type	Test weights
Product group	Safety sets
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761245716
Planned start date	2019-1HJ
Customs tariff number	84239090

Optional Services

DAkkS Calibration	962-400-600
-------------------	-------------

Technical data - Measuring system

Nominal weights, included	100 g (M1); 1000 g (F1)
Nominal value	100 g & 1 kg
OIML class	F1;M1

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	255×135×140 mm
Gross weight	3 kg
Net weight	1.2 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	1.8 kg
Delivery time	1 d

Technical data - Verification

Verification possible	Yes
-----------------------	-----

Functions

KERN 315-211-407

Safety Set 100 g (M1); 1000 g (F1)



Standard

