

KERN 366-96

Individual weight, OIML Class M3, 5 kg, Knob, cast iron lacquered



Article codes, product groups & Status information

Model article-code	366-96
Model-Serie	366
Assortment	KERN
Product Type	Test weights
Product group	OIML M3 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761482104
Predecessor Model	366-96-2.024a
Planned start date	2025-2HJ
Customs tariff number	84239090

Technical data - Primary

Nominal value	5 kg
Tolerance (OIML)	2,5 g

Technical data - Construction

Material	Cast iron lacquered
----------	---------------------

Technical data - Measuring system

Nominal value	5 kg
OIML class	M3
OIML conform	Yes
Tolerance (OIML)	2,5 g
Density (kg/m ³)	7.100
Volume (cm ³)	704,2254

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	94×94×146 mm
Gross weight	5,0 kg
Net weight	5,0 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	No
-----------------------	----

Optional Services

DAkKS Calibration	962-643
-------------------	---------

Functions

KERN 366-96

Individual weight, OIML Class M3, 5 kg, Knob, cast iron lacquered



Standard



Option

