

PRODUCT BROCHURE

KERN 316-12

Date: 19.09.2025

ENGLISH



Individual weight, OIML Class E2, 2 kg, Compact form, stainless steel polished

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

- Individual weight, OIML Class E2, 2 kg, Compact form, stainless steel polished (OIML)
- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	316-12
Model-Serie	316
Assortment	KERN
Product Type	Test weights
Product group	OIML E2 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761094932
Customs tariff number	84239090

Technical data - Primary

Nominal value	2 kg
Tolerance (OIML)	3 mg
Case aluminium	317-120-600
Case plastic	317-120-400
Case wood	317-120-100

Technical data - Construction

Material	stainless steel polished
----------	--------------------------

Technical data - Measuring system

Nominal value	2 kg
OIML class	E2
OIML conform	yes
Tolerance (OIML)	3 mg
Density (kg/m ³)	7950
Volume (cm ³)	250

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	60×60×94 mm
Gross weight	2 kg
Net weight	2 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	no
-----------------------	----

Optional Services

DAkkS Calibration	962-342
-------------------	---------

FUNCTIONS

Standard



ACCESSORIES

Model	Description
317-280	Gloves KERN 317-280
317-290	Gloves KERN 317-290
317-120-100	Wooden Box for Set of Weights KERN 317-120-100
317-120-400	Plastic Weight Box KERN 317-120-400
317-120-600	Aluminium Protective Box KERN 317-120-600

SERVICES

Model	Description
962-342	DAkkS Calibration KERN 962-342
962-342R	DAkkS Recalibration KERN 962-342R
921-342	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-342
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351

RELATED ARTICLES DATA

Model name	Nominal value	Tolerance (OIML)	Case aluminium	Case plastic	Case wood
316-01	1 g	0,03 mg	317-010-600	317-020-400	317-010-100
316-02	2 g	0,04 mg	317-020-600	317-020-400	317-020-100
316-03	5 g	0,05 mg	317-030-600	317-030-400	317-030-100
316-04	10 g	0,06 mg	317-040-600	317-040-400	317-040-100
316-05	20 g	0,08 mg	317-050-600	317-050-400	317-050-100
316-06	50 g	0,1 mg	317-060-600	317-060-400	317-060-100
316-07	100 g	0,16 mg	317-070-600	317-070-400	317-070-100
316-08	200 g	0,3 mg	317-080-600	317-080-400	317-080-100
316-09	500 g	0,8 mg	317-090-600	317-090-400	317-090-100
316-11	1 kg	1,6 mg	317-110-600	317-110-400	317-110-100
316-12	2 kg	3 mg	317-120-600	317-120-400	317-120-100

Model name	Nominal value	Tolerance (OIML)	Case aluminium	Case plastic	Case wood
316-13	5 kg	8 mg	317-130-600	317-130-400	317-130-100
316-14	10 kg	16 mg	317-140-600	317-140-400	317-140-100