

PRODUCT BROCHURE

KERN 318-226

Date: 10.03.2026

ENGLISH



Set of milligram weights, E2, 1 mg - 500 mg, Platelet, stainless steel polished, in aluminium case

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

PRODUCT SHOWCASE



SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	318-226
Model-Serie	318
Assortment	KERN
Product Type	Test weights
Product group	OIML E2 Test weights
Product usage type	Main article
EAN code (Model)	4045761223615
Customs tariff number	84239090
Country of Origin	CN

Technical data - Primary

Nominal value	1 mg - 500 mg
Case aluminium	313-010-600
Weight set case	in aluminium case

Technical data - Construction

Material	stainless steel polished
Is weight set	yes
Number of weights in set	13
Type of case	in aluminium case

Technical data - Measuring system

Nominal value	1 mg - 500 mg
OIML class	E2
OIML conform	yes

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	165×158×37 mm
Gross weight	0.40 kg
Net weight	0.35 kg
Shipping method	Parcel service

Packaging component - 80 g
by weight - plastics

Delivery time 1 d

Optional Services

DAkkS Calibration 962-350

FUNCTIONS

Standard



Option



SERVICES

Model	Description
962-350	DAkkS Calibration KERN 962-350
962-350R	DAkkS Recalibration KERN 962-350R
921-350	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-350
920-701	Marking (Last digits of serial number) KERN 920-701

SPARE PARTS

Model	Description
313-010-600	Aluminium Protective Box KERN 313-010-600
315-242	Tweezers KERN 315-242
318-270	Dusting Brush KERN 318-270